Blanchard Valley Health System Location:	_
Training / Competency As sessment for Laboratory Point of Care Testing	Date of Training:
Method: POCT Medline 10SG Urinalysis Reagent Strips on Medline 120 Mi	ini Urine Analyzer
Test analytes: Urine Qualitative Dipstick Leukocytes, Nitrite, Urobilinogen	n, Protein, pH, Blood, SG, Ketone, Bilirubin, Glu
Associate Name and Title:	_
(print legibly) Associate ID number: Supervisor:	
CLIA Complexity: Waived (2 areas must be assessed) (CAP approved areas	1 – 6) (JC approved areas 1, 3, 5, 6)
Training assessment after initial training (new hires) and/or new method and	nd/or new test analytes Date:
Competency assessment annually	Date Due:
123 Perform QC testing as required. (Properly running acceptable test results meets CAP/JC competency requirements for elements both levels of liquid urine controls on analyzer. Attach printout. Results ac	cceptable? Yes / No PASS / FAIL (circle)
4 Perform daily cleaning. (Properly performing instrument mainten observation and documenting acceptable results meets CAP/JC competer	
Clean strip holder and the White Calibration Circle as necessary using a cotton	swab/ball with distilled water.
Dry with a clean, dry cotton sw ab/ball.	PASS / FAIL (cirde)
5 Perform testing on an unknown sample as required. (Analyzing a results meets CAP/JC competency requirements for element 5.)	
UNKNOWN Test: (Patient or control or CAP Proficiency sample) Test sample	on a nalyzer. Attach printout.
COLORCLARITY / APPEARANCE	PASS / FAIL (circle)
Associate Signature:	Date:
Signature of Qualified Trainer/Assessor:	Date:
Job title of Qualified Trainer/Assessor:	Satisfactory / Unsatisfactory
Reviewed by:	

Date: