

BLANCHARD VALLEY PHLEBOTOMY SPECIMEN REQUIREMENTS

HEP GEL

(LIGHT GREEN/GEL)

ALBUMIN
 ACETAMINOPHEN
 ALK PHOS
 ALT (SGPT)
 AMYLASE
 AST (SGOT)
 B12
 BUN
 CALCIUM
 CHLORIDE
 *CHOLESTROL
 CO2
 CORTISOL
 CPK
 CREATININE
 CRP
 DIGOXIN**Findlay**
 DIRECT BILI
 ESTRADIOL
 FERRITIN
 FOLATE (FOLIC ACID)
 FREE T3 AND FREE T4
 FSH
 GAM (IgG/IgA/IgM)
 GENTAMICIN
 GGT
 *GLUCOSE
 *HDL
 IRON
 LDH
 LH
 LIPASE
 MAGNESIUM
 MYOGLOBIN
 PHOSPHORUS
 POTASSIUM
 PREALBUMIN
 PROLACTIN
 SALICYLATE
 SODIUM
 TOTAL T3
 T3 UPTAKE
 T4
 TOTAL
 TESTOSTERONE
 TIBC (TOTAL IRON)
 TOBRAMYCIN
 TOTAL BILI
 TOTAL PROTEIN
 TRANSFERRIN

*TRIGLYCERIDE
 TROPONIN
 TSH
 URIC ACID
 VANCOMYCIN
 VITAMIN D-25
 HYDROXY

CHEMISTRY

PROFILES: (LIGHT GREEN/GEL)

*BASIC (BMP):

SODIUM
 POTASSIUM
 CHLORIDE
 CO2
 GLUCOSE
 BUN
 CREATININE
 CALCIUM

*COMPREHENSIVE

(CMP):

SODIUM
 POTASSIUM
 CHLORIDE
 CO2
 GLUCOSE
 BUN
 CREATININE
 CALCIUM
 ALK PHOS
 AST (SGPT)
 ALT (SGOT)
 TOTAL PROTEIN
 ALBUMIN
 TOTAL BILI

HEPATIC PANEL:

ALK PHOS
 ALT (SGPT)
 AST (SGOT)
 TOTAL PROTEIN
 TOTAL BILI
 DIRECT BILI
 INDIRECT BILI
 ALBUMIN

RENAL PANEL :

SODIUM
 POTASSIUM

CHLORIDE
 CO2
 GLUCOSE
 BUN
 CREATININE
 CALCIUM
 Albumin
 Phosphorus

*LIPID PROFILE:

TRIGLYCERIDE
 TOTAL CHOLESTROL
 HDL
 LDL
 CARDIAC RISK
 FACTOR

PINK (EDTA)

ABORH
 ANTIBODY SCREEN
 CROSSMATCH

PLAIN RED NO GEL

BETA HCG
 BHOB
 DILANTIN
 DIGOXIN**Bluffton
 PHENOBARBITAL
 PROGESTERONE
 SERUM PREGNANCY
 TEGRETOL
 VALPROIC ACID

Yellow (SST)

CA125
 CA 15-3
 CEA
 ELECTROPHORESIS
 HEPATITIS
 IMMUNOFIXATION
 LITHIUM
 OSMO
 PSA
 FREE PSA

LAVENDER (EDTA)

BNP
 CBC
 ESR (SED RATE)
 HBA1C
 HIV
 HHP

INSULIN

BLUE (CITRATE)

DIMER
 PTT
 PTINR
 FIBRINOGEN
 HYPERCOAG PANEL:
 4 BLUE CITRATE
 TUBES
 1 HEP GEL (ON ICE)
 1 ACD TUBE

GREY

ALCOHOL (ETOH)
 Lactic Acid

SPECIAL PREP

AMMONIA:

- PLAIN GREEN (LI HEP)
- ON ICE
- SPUN 15 MINS

LACTIC ACID:

- GREY
- ON ICE
- SPUN 15 MINS

TESTOSTERONE

FREE & TOTAL:

- PLAIN RED

IONIZED CALCIUM:

- SST
- TUBE MUST BE FULL
- DO NOT UNCAP

RENIN ACTIVITY:

- PRECHILLED EDTA
- SEND TO LAB ON ICE

(*) REQUIRES FASTING