

## BLANCHARD VALLEY PHLEBOTOMY SPECIMEN REQUIREMENTS

### HEP GEL

#### (LIGHT GREEN/GEL)

ALBUMIN  
 ACETAMINOPHEN  
 ALK PHOS  
 ALT (SGPT)  
 AMYLASE  
 AST (SGOT)  
 B12  
 BUN  
 CALCIUM  
 CHLORIDE  
 \*CHOLESTROL  
 CO2  
 CORTISOL  
 CPK  
 CREATININE  
 CRP  
 DIGOXIN\*\*Findlay\*\*  
 DIRECT BILI  
 ESTRADIOL  
 FERRITIN  
 FOLATE (FOLIC ACID)  
 FREE T3 AND FREE T4  
 FSH  
 GAM (IgG/IgA/IgM)  
 GENTAMICIN  
 GGT  
 \*GLUCOSE  
 \*HDL  
 IRON  
 LDH  
 LH  
 LIPASE  
 MAGNESIUM  
 MYOGLOBIN  
 PHOSPHORUS  
 POTASSIUM  
 PREALBUMIN  
 PROLACTIN  
 SALICYLATE  
 SODIUM  
 TOTAL T3  
 T3 UPTAKE  
 T4  
 TOTAL  
 TESTOSTERONE  
 TIBC (TOTAL IRON)  
 TOBRAMYCIN  
 TOTAL BILI  
 TOTAL PROTEIN  
 TRANSFERRIN

\*TRIGLYCERIDE  
 TROPONIN  
 TSH  
 URIC ACID  
 VANCOMYCIN  
 VITAMIN D-25  
 HYDROXY

### CHEMISTRY

#### PROFILES: (LIGHT GREEN/GEL)

#### \*BASIC (BMP):

SODIUM  
 POTASSIUM  
 CHLORIDE  
 CO2  
 GLUCOSE  
 BUN  
 CREATININE  
 CALCIUM

#### \*COMPREHENSIVE

#### (CMP):

SODIUM  
 POTASSIUM  
 CHLORIDE  
 CO2  
 GLUCOSE  
 BUN  
 CREATININE  
 CALCIUM  
 ALK PHOS  
 AST (SGPT)  
 ALT (SGOT)  
 TOTAL PROTEIN  
 ALBUMIN  
 TOTAL BILI

#### HEPATIC PANEL:

ALK PHOS  
 ALT (SGPT)  
 AST (SGOT)  
 TOTAL PROTEIN  
 TOTAL BILI  
 DIRECT BILI  
 INDIRECT BILI  
 ALBUMIN

#### RENAL PANEL :

SODIUM  
 POTASSIUM

CHLORIDE  
 CO2  
 GLUCOSE  
 BUN  
 CREATININE  
 CALCIUM  
 Albumin  
 Phosphorus

#### \*LIPID PROFILE:

TRIGLYCERIDE  
 TOTAL CHOLESTROL  
 HDL  
 LDL  
 CARDIAC RISK  
 FACTOR

#### PINK (EDTA)

ABORH  
 ANTIBODY SCREEN  
 CROSSMATCH

#### PLAIN RED NO GEL

BETA HCG  
 BHOB  
 DILANTIN  
 DIGOXIN\*\*Bluffton  
 PHENOBARBITAL  
 PROGESTERONE  
 SERUM PREGNANCY  
 TEGRETOL  
 VALPROIC ACID

#### Yellow (SST)

CA125  
 CA 15-3  
 CEA  
 ELECTROPHORESIS  
 HEPATITIS  
 IMMUNOFIXATION  
 LITHIUM  
 OSMO  
 PSA  
 FREE PSA

#### LAVENDER (EDTA)

BNP  
 CBC  
 ESR (SED RATE)  
 HBA1C  
 HIV  
 HHP

INSULIN

#### BLUE (CITRATE)

DIMER  
 PTT  
 PTINR  
 FIBRINOGEN  
 HYPERCOAG PANEL:  
 4 BLUE CITRATE  
 TUBES  
 1 HEP GEL (ON ICE)  
 1 ACD TUBE

#### GREY

ALCOHOL (ETOH)  
 Lactic Acid

#### SPECIAL PREP

#### AMMONIA:

- PLAIN GREEN (LI HEP)
- ON ICE
- SPUN 15 MINS

#### LACTIC ACID:

- GREY
- ON ICE
- SPUN 15 MINS

#### TESTOSTERONE

#### FREE & TOTAL:

- PLAIN RED

#### IONIZED CALCIUM:

- SST
- TUBE MUST BE FULL
- DO NOT UNCAP

#### RENIN ACTIVITY:

- PRECHILLED EDTA
- SEND TO LAB ON ICE

(\*) REQUIRES FASTING