








SWAB REFERENCE GUIDE

Specimen Type	General Tests Associated	Photo
<p><b>E Swab</b></p>	<ul style="list-style-type: none"> <li>• Cultures, aerobic and anaerobic</li> <li>• MRSA</li> <li>• Group B Strep</li> <li>• Throat Culture</li> </ul>	
<p><b>Nasal Swab with Universal Transport Media (M4)</b></p>	<ul style="list-style-type: none"> <li>• Influenza A/B (PCR)</li> <li>• RSV (PCR)</li> <li>• SARS-CoV-2 RNA by RT-PCR</li> </ul>	
<p><b>Sterile Foam Tip (Dry)</b> Covid Swabs must be tested within 1 hour **Do NOT place swab back in original packaging, use plastic screw cap tube for transport.</p>	<ul style="list-style-type: none"> <li>• SARS-COVID19 RNA Rapid (nasal)</li> <li>• SARS-CoV-2 Antigen (nasal)</li> <li>• Influenza A/B RNA Rapid (nasal)</li> <li>• Strep A DNA Rapid (throat)</li> </ul>	
<p><b>Nasopharyngeal Swab</b> Flu and RSV swabs must be tested within 4 hours</p>	<ul style="list-style-type: none"> <li>• Influenza A/B RNA Rapid</li> <li>• RSV RNA Rapid</li> </ul>	

SWAB REFERENCE GUIDE

<p><b>Aptima® Multitest Swab</b> Specimen Collection Kit for Vaginal, Throat, and Rectal Swab Specimens</p>	<ul style="list-style-type: none"> <li>• GC</li> <li>• Chlamydia</li> </ul>	
<p><b>Aptima® Unisex Swab Specimen</b> Collection Kit for Endocervical and Male Urethral Swab Specimens</p>	<ul style="list-style-type: none"> <li>• GC</li> <li>• Chlamydia</li> </ul>	
<p><b>Red Duo Swab</b></p>	<ul style="list-style-type: none"> <li>• Culture Backup Strep ***Bluffton ONLY***</li> </ul>	
<p><b>Xpert Swab Specimen Collection Kit</b></p>	<ul style="list-style-type: none"> <li>• Vaginal Panel DNA (MVP)</li> </ul>	