





SWAB REFERENCE GUIDE

Specimen Type	General Tests Associated	Photo
<p>E Swab</p>	<ul style="list-style-type: none"> • Cultures, aerobic and anaerobic • MRSA • Group B Strep • Throat Culture 	
<p>Nasal Swab with UTM</p>	<ul style="list-style-type: none"> • Influenza A/B (PCR) • RSV (PCR) • SARS-CoV-2 RNA by RT-PCR 	
<p>Sterile Foam Tip (Dry) Covid Swabs must be tested within 1 hour **Do NOT place swab back in original packaging, use plastic screw cap tube for transport.</p>	<ul style="list-style-type: none"> • SARS-COVID19 RNA Rapid (nasal) • SARS-CoV-2 Antigen (nasal) • Influenza A/B RNA Rapid (nasal) • Strep A DNA Rapid (throat) 	
<p>Nasopharyngeal Swab Flu and RSV swabs must be tested within 4 hours</p>	<ul style="list-style-type: none"> • Influenza A/B RNA Rapid • RSV RNA Rapid 	

SWAB REFERENCE GUIDE

<p>Red Duo Swab</p>	<ul style="list-style-type: none"> • Culture Backup Strep ***Bluffton ONLY*** 	
<p>Affirm Swab</p>	<ul style="list-style-type: none"> • Affirm DNA Panel 	
<p>Dry Culture Swab **See kit for specific instructions</p>	<ul style="list-style-type: none"> • Monkeypox 	