## ID NOW™ COVID-19 2.0 External Quality Control Log

External QC testing is recommended:

- When a new shipment of kits is received
- When a new untrained operator performs testing
- The first time an assay is run on an instrument
- When required by local, state, and/or federal regulations, accrediting groups, or your lab's Quality Control procedures

ID NOW INSTRUMENT SERIAL NUMBER



DATE	ID NOW™ COVID-19 2.0 KIT LOT/EXP	POSITIVE CONTROL LOT/EXP	NEGATIVE CONTROL LOT/EXP	POSITIVE CONTROL RESULT	NEGATIVE CONTROL RESULT

Supervisory Review by:	Date:

FOR EXTERNAL USE. PRINT AND DISTRIBUTION PERMITTED.

© 2022 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. 120004477-05 Rev.02 5/2022



