



**ID NOW™ Influenza A & B 2 External Quality Control**

Form revised 8/9/2023 LTR

External QC testing is recommended:

- When a new shipment of kits is received
- When a new untrained operator performs testing
- The first time an assay is run on an instrument
- When required by local, state, and/or federal regulations, accrediting groups, or your lab's Quality Control procedures.

ID Now INSTRUMENT SERIAL NUMBER \_\_\_\_\_

DATE	ID NOW™ Influenza A & B 2 KIT LOT/EXP	POSITIVE CONTROL LOT/EXP	NEGATIVE CONTROL LOT/EXP	POSITIVE CONTROL RESULT	NEGATIVE CONTROL RESULT	TESTER'S INITIALS	COMMENTS
R E							

Supervisory Review by: \_\_\_\_\_

Date: \_\_\_\_\_

FOR EXTERNAL USE. PRINT AND DISTRIBUTION PERMITTED.

© 2022 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

120004477-05 Rev.02 5/2022