



ID NOW™ Quality Control Printer Labels

External positive and negative controls, which are included in the kit, should be tested when an assay is run on the instrument for the first time, once with each new shipment and once for each untrained operator, following a software upgrade, or in order to conform to state, local and/or federal regulations, accrediting groups, or your laboratory's standard Quality Control procedures.

KIT REC'D DATE:	KIT EXP. DATE	POS QC EXP. DATE	NEG QC EXP. DATE	TECH	DATE
POSITIVE CONTROL RESULT LABEL		NEGATIVE CONTROL RESULT LABEL			
<i>Please initial QC label if applicable</i>		<i>Please initial QC label if applicable</i>			
KIT REC'D DATE:	KIT EXP. DATE	POS QC EXP. DATE	NEG QC EXP. DATE	TECH	DATE
POSITIVE CONTROL RESULT LABEL		NEGATIVE CONTROL RESULT LABEL			
<i>Please initial QC label if applicable</i>		<i>Please initial QC label if applicable</i>			
COMMENTS:					