

Method: **POCT Nova StatStrip Xpress 2**

QC and Maintenance

Year:

Test analyte: **Whole Blood Glucose**

Month:

Dept. / Location:

Meter SN:

Level 1 Lot #		Exp Date		Range					
Level 3 Lot #		Exp Date		Range					
Test Strips Lot#		Exp Date		Mthly batteries check Acceptable? YES / NO					
				Date & Initials:					
DAY	GLUCOSE METER NOT IN USE (Initials)	OPEN EXP DATE TEST STRIPS (6 mos) & CONTROLS (3 mos) YES/NO (Circle)		BATTERIES REPLACED YES/NO (Circle)	Level 1 LOW CONTROL result	Level 3 HIGH CONTROL result	Acceptable? YES/NO (Circle)	Corrective action taken	Trained operator initials
1		Y	N	Y	N		Y	N	
2		Y	N	Y	N		Y	N	
3		Y	N	Y	N		Y	N	
4		Y	N	Y	N		Y	N	
5		Y	N	Y	N		Y	N	
6		Y	N	Y	N		Y	N	
7		Y	N	Y	N		Y	N	
8		Y	N	Y	N		Y	N	
9		Y	N	Y	N		Y	N	
10		Y	N	Y	N		Y	N	
11		Y	N	Y	N		Y	N	
12		Y	N	Y	N		Y	N	
13		Y	N	Y	N		Y	N	
14		Y	N	Y	N		Y	N	
15		Y	N	Y	N		Y	N	
16		Y	N	Y	N		Y	N	
17		Y	N	Y	N		Y	N	
18		Y	N	Y	N		Y	N	
19		Y	N	Y	N		Y	N	
20		Y	N	Y	N		Y	N	
21		Y	N	Y	N		Y	N	
22		Y	N	Y	N		Y	N	
23		Y	N	Y	N		Y	N	
24		Y	N	Y	N		Y	N	
25		Y	N	Y	N		Y	N	
26		Y	N	Y	N		Y	N	
27		Y	N	Y	N		Y	N	
28		Y	N	Y	N		Y	N	
29		Y	N	Y	N		Y	N	
30		Y	N	Y	N		Y	N	
31		Y	N	Y	N		Y	N	

Note: CAP/JC sites must enter results in Cerner too.

Monthly supervisory review by: _____

Date: _____